INSTRUCTION FOR USE

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CORITON Medical Device Support Systems
医疗推车、移动支架、支臂及上墙支架

CORITON_IFU_2026 | Revision 03 | Issue 2026-07

English (EN) — 2	العربية (AR) — 3	Български (BG) — 4	Čeština (CS) — 5	Dansk (DA) — 6
Deutsch (DE) — 7	Ελληνικά (EL) — 8	Español (ES) — 9	Eesti (ET) — 10	Suomi (FI) — 11
Français (FR) — 12	Hrvatski (HR) — 13	Magyar (HU) — 14	Italiano (IT) — 15	日本語 (JA) — 16
Lietuvių (LT) — 17	Latviešu (LV) — 18	Nederlands (NL) — 19	Norsk (NO) — 20	Polski (PL) — 21
Português (PT) — 22	Română (RO) — 23	Русский (RU) — 24	Slovenčina (SK) — 25	Slovenščina (SL) — 26
Srpski (SR) — 27	Svenska (SV) — 28	Türkçe (TR) — 29	Українська (UA) — 30	中文 (ZH) — 31

INTENDED USE

CORITON medical carts and mounting arms are designed to hold, position, and transport monitors, medical devices, keyboards, laptops, and other medical equipment (hereinafter "ME equipment"). The systems can be configured with various accessories including shelves, drawers, monitor mounts, keyboard trays, and cable management solutions. CORITON systems support medical applications, but do not serve an independent medical purpose on their own.

This instruction for use is intended for the operator and user of CORITON systems. Instructions for ME equipment and their accessories are not included herein.

PRODUCT IDENTIFICATION AND TERMINOLOGY

Depending on the model, configuration and market, the product label may use Cart, Trolley, Roll Stand, Arm, Mounting Arm, Wall Mount or Wall-mounted Arm. This IFU collectively refers to them as CORITON support systems. Always follow the label, model/UDI and model-specific installation, load, compatibility and use limits.

SAFETY INSTRUCTIONS

- Do not lean on CORITON carts or mounting arms. This may damage equipment or cause instability. Risk of injury.
- Avoid ingress of disinfectants and cleaning liquids into the system interior.
- Do not use steel wool or abrasive agents.
- Compatibility of disinfectants or cleaning fluids should be tested beforehand on a non-critical area.
- Wipe external surfaces with a soft, lint-free cloth and a healthcare surface cleaner or disinfectant compatible with stainless steel, anodized aluminium, powder-coated surfaces and ABS. Follow the agent manufacturer's instructions.
- Use of acetone and trichloroethylene will cause permanent damage to surfaces.
- Cables integrated into the CORITON system may also be affected by unsuitable agents. In case of doubt, contact the manufacturer.
- Never pull on cables or wires. Damaged cables, lines, or defective ME equipment must be immediately disconnected from the power supply.
- The rotation range is limited by internal rotation stops. Do not attempt to force beyond the rotation stop, as this may cause damage to cables or equipment and endanger operators.
- The rotation range of CORITON systems must be outside the patient area. No other devices may be located within the arm's rotation range. Risk of damage!
- Move carefully, and mind pinch points.
- During transport, secure all support arms in a safe position.
- Always move carts with both hands to prevent tipping.
- Installed ME equipment or cables may operate at hazardous electrical voltages.
- Do not connect live cables to the power supply until installation is complete. Before performing inspection/disassembly, disconnect all live cables from the power supply and equipment.
- Before use, read and follow this instruction for use, including all safety information and operating instructions.

MAINTENANCE

- CORITON systems are maintenance-free.
- Inspection should be performed by qualified personnel at least once a year.

LOCKING MECHANISMS

Where applicable, CORITON systems are equipped with locking levers/parking brakes. To change the position, secure the ME equipment by hand, release the locking lever, re-adjust the position, and re-tighten the locking lever/parking brake.

Carts are equipped with one or more lockable casters. After positioning, lock all brakes by engaging the lever with your foot.

ASSEMBLY / INSTALLATION

Assembly, disassembly, cabling and integration of ME equipment, as well as inspection, may only be performed by qualified personnel in accordance with the corresponding assembly instructions.

Ensure that the selected fastening option or mounting adapter can withstand the respective weight load. Use of suitable cables and implementation of required grounding measures are the responsibility of the customer/system integrator.

CORITON recommends consulting a safety officer and biomedical technician regarding installation and use of CORITON systems.

CORITON mounting instructions with inspection plan and instructions for use are available at:

<https://www.coriton.com/services/media-center/>

PARTS AND ADAPTERS: Purely mechanical, interface-conversion or cable-management parts are not medical devices/accessories under Regulation (EU) 2017/745; parts supporting a medical function may be accessories. Follow the label and compatibility, load and installation limits.

CE MARKING

CORITON systems and their accessories comply with applicable safety requirements.

SERIOUS INCIDENT REPORTING


All serious incidents occurring in connection with or caused by CORITON systems must be reported to the manufacturer and the competent authorities.

MANUFACTURER

CORITON Instrument Co., Ltd
602, Bldg 1, No.1 Danqing Rd, Innovation Mfg Park,
Kengzi, Pingshan, Shenzhen 518122, China
info@coriton.com www.coriton.com

EU AUTHORISED REPRESENTATIVE

Umedwings Netherlands B.V.
Treubstraat 1, 2288EG, Rijswijk, The Netherlands
ar@umedwings.eu

UK RESPONSIBLE PERSON

Umedwings UK Ltd
The Foundry, 2 Smiths Square, 77 Fulham Palace Road,
London, W6 8AF
ukrp@umedwings.uk

RECYCLING INFORMATION

CORITON systems contain valuable raw materials and can be largely recycled. Electrical equipment must be collected and disposed of separately.

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